



YOUR LAB NAME

123 Main Street

Anywhere, NC USA 28273

CLIA ID : 34D9999999 Lab Director : Dr. Sample, PhD

PATIENT NAME	DOB	GENDER	ACCESSION#	COLLECTED	RECEIVED
Test T Test	08/05/1981	Male	CV20-1049	12/08/2020	12/08/2020

ORDERING PROVIDER	LOCATION	REPORTED
COVID19 Provider	ABC Diagnostics	12/08/2020

**Test Summary**

**COVID-19 (Saliva)** **Notes:**  
**Sample Type:** Saliva

**Result Summary: Tests Performed, Non-Detected**

Test Name	Result	Comments
<b>COVID-19 (Saliva)</b>		
SARS-CoV-2	INCONCLUSIVE	Unable to determine result. Retest patient if clinically indicated.

**Methodologies** : Saliva samples are pretreated with Proteinase K prior to Reverse Transcriptase Polymerase Chain Reaction (RT-PCR) using the CDC primer probes (CDC 2019-nCoV Real-Time RT-PCR Diagnostic Panel), Taqpath one step master mix and Twist synthetic RNA control.

**Limitations** : False-negative results may arise from improper sample collection, degradation of the SARS-CoV-2 RNA during shipping/storage, specimen collection after SARS-CoV-2 RNA can no longer be found in the specimen matrix, using unauthorized assay reagents, the presence of RT-PCR inhibitors, mutation in the SARS-CoV-2 virus, failure to follow instructions for use. False-positive results may arise from cross contamination during specimen handling or preparation, cross contamination between patient samples, specimen mix-up – RNA contamination during product handling.

SAMPLE

**Disclaimer:** This test has been developed and its performance characteristics determined by Your Lab Name . It has not been cleared or approved by the FDA. The laboratory is regulated under CLIA as qualified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.



YOUR LAB NAME

123 Main Street

Anywhere USA 28273

CLIA ID : 34D9999999 Lab Director : Dr. Sample, PhD

PATIENT NAME	DOB	GENDER	ACCESSION#	COLLECTED	RECEIVED
Test T Test	08/05/1981	Male	CV20-1049	12/08/2020	12/08/2020

ORDERING PROVIDER	LOCATION	REPORTED
COVID19 Provider	ABC Diagnostics	01/19/2021

**Test Summary**

**COVID-19 (Saliva)**

**Notes:**

**Sample Type:** Saliva

**Result Summary: Tests Performed, Non-Detected**

Test Name	Result	Comments
COVID-19 (Saliva)		
SARS-CoV-2	REJECTED	

**Methodologies** : Saliva samples are pretreated with Proteinase K prior to Reverse Transcriptase Polymerase Chain Reaction (RT-PCR) using the CDC primer probes (CDC 2019-nCoV Real-Time RT-PCR Diagnostic Panel), Taqpath one step master mix and Twist synthetic RNA control.

**Limitations** : False-negative results may arise from improper sample collection, degradation of the SARS-CoV-2 RNA during shipping/storage, specimen collection after SARS-CoV-2 RNA can no longer be found in the specimen matrix, using unauthorized assay reagents, the presence of RT-PCR inhibitors, mutation in the SARS-CoV-2 virus, failure to follow instructions for use. False-positive results may arise from cross contamination during specimen handling or preparation, cross contamination between patient samples, specimen mix-up – RNA contamination during product handling.

SAMPLE

**Disclaimer:** This test has been developed and its performance characteristics determined by YOUR LAB NAME . It has not been cleared or approved by the FDA. The laboratory is regulated under CLIA as qualified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.



YOUR LAB NAME

123 Main Street

Anywhere, NC USA 28273

CLIA ID : 34D9999999 Lab Director : Dr. Sample, PhD

PATIENT NAME	DOB	GENDER	ACCESSION#	COLLECTED	RECEIVED
Test T Test	08/05/1981	Male	CV20-1048	12/07/2020	12/07/2020

ORDERING PROVIDER	LOCATION	REPORTED
Test Provider	ABC Diagnostics	12/08/2020

**Test Summary**

**SARS-CoV-2**

**Notes:**

**Sample Type:** Saliva

**Result Summary: Tests Performed, Non-Detected**

Test Name	Result	Comments
<b>SARS-CoV-2</b>		
SARS-CoV-2	NOT DETECTED	

**Methodologies :** Saliva samples are pretreated with Proteinase K prior to Reverse Transcriptase Polymerase Chain Reaction (RT-PCR) using the CDC primer probes (CDC 2019-nCoV Real-Time RT-PCR Diagnostic Panel), Taqpath one step master mix and Twist synthetic RNA control.

**Limitations :** False-negative results may arise from improper sample collection, degradation of the SARS-CoV-2 RNA during shipping/storage, specimen collection after SARS-CoV-2 RNA can no longer be found in the specimen matrix, using unauthorized assay reagents, the presence of RT-PCR inhibitors, mutation in the SARS-CoV-2 virus, failure to follow instructions for use. False-positive results may arise from cross contamination during specimen handling or preparation, cross contamination between patient samples, specimen mix-up – RNA contamination during product handling.

SAMPLE

**Disclaimer:** This test has been developed and its performance characteristics determined by YOUR LAB NAME . It has not been cleared or approved by the FDA. The laboratory is regulated under CLIA as qualified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.



YOUR LAB NAME

9123 Main Street

Anywhere, NC USA 28273

CLIA ID : 34D9999999 Lab Director : Dr. Sample, PhD

PATIENT NAME	DOB	GENDER	ACCESSION#	COLLECTED	RECEIVED
Test T Test	08/05/1981	Male	CV20-1190	12/28/2020	12/29/2020

ORDERING PROVIDER	LOCATION	REPORTED
Test Provider	ABC Diagnostics	12/29/2020

**Test Summary**

**COVID-19 (Saliva)**

**Notes:**

**Sample Type:** Saliva

**Result Summary: Tests Performed, Detected**

Test Name	Result	Comments
<b>COVID-19 (Saliva)</b>		
SARS-CoV-2	DETECTED	NOTIFICATION REQUIRED

**Methodologies :** Saliva samples are pretreated with Proteinase K prior to Reverse Transcriptase Polymerase Chain Reaction (RT-PCR) using the CDC primer probes (CDC 2019-nCoV Real-Time RT-PCR Diagnostic Panel), Taqpath one step master mix and Twist synthetic RNA control.

**Limitations :** False-negative results may arise from improper sample collection, degradation of the SARS-CoV-2 RNA during shipping/storage, specimen collection after SARS-CoV-2 RNA can no longer be found in the specimen matrix, using unauthorized assay reagents, the presence of RT-PCR inhibitors, mutation in the SARS-CoV-2 virus, failure to follow instructions for use. False-positive results may arise from cross contamination during specimen handling or preparation, cross contamination between patient samples, specimen mix-up – RNA contamination during product handling.

SAMPLE

**Disclaimer:** This test has been developed and its performance characteristics determined by YOUR LAB NAME . It has not been cleared or approved by the FDA. The laboratory is regulated under CLIA as qualified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.